

PioNext Mini

All-in-One Dental 3D Printing Solution



Miniaturized
Design



10min
Fast Printing



Integrated
Printing & Curing



2 times
Screen Lifespan



Wash Mini



ACF Fast Vat



HDF Clear Vat



One machine handles the entire production process for
true chairside efficiency

35-Minute Delivery

Ceramic Crown Chairside Workflow



Learn More

One machine handles the entire production process for true chairside efficiency



More Applications



Technical Specifications

Printing Size	80mm×70mm×80mm	Resolution	9K (8520*4320px)
Machine Size	200mm×170mm×241mm	Pixel Size	14*19μm
Package Size	310mm×307mm×340mm	Layer Thickness	0.01-0.1mm
Curing Size	76.5mm×72.5mm×37mm	Touch Screen	3.97-inch Capacitive Color Touch Screen
Light Intensity(Curing)	> 150 mW/cm ² , 405nm	Printing Speed	Max.60mm/h
N.W.	4.86kg	Slicing Software	PioNext BOX / CHITUBOX Dental
G.W.	6.8kg	Data Transmission	USB/WIFI
Printing Technology	LCD	Input Voltage	100-240V~, 50/60Hz
Light Source Wavelength	385nm、405nm	Power	120W
Supported Consumables	Ceramic Crown Resin / Surgical Guide Resin / Crown & Bridge Resin / Temporary Restoration Resin / Denture Base Resin / Removable Die Model Resin / Ortho Model Resin / Ortho Model 2.0 Resin and other official authorized resin		

PioNext Mini

All-in-One Dental 3D Printing Solution

