


White ENDOSEAL Bioceramic Sealer

Beyond Sealing & Filling

Not just a sealer! Not just a filler!



What Makes White Endoseal Bioceramic Sealer Stand Out?


In the world of bioceramic sealers, bioactivity is essential, but achieving a truly stable apical seal demands more.

White Endoseal Bioceramic Sealer sets apart from the rest because it sets quickly, has a high washout resistance and is more radiopaque compared to current market leaders.

A combination of three patented effective ingredients that regulate immune response, promote cell differentiation and enhance healing efficiency.

Tricalcium silicate, Lithium carbonate, Dimethyl sulfoxide

White Endoseal Bioceramic Sealer presents a new treatment method for root canal treatment with healing and regeneration, rather than the conventional treatment method based on repair.



Key Features

White ROOT CANAL FILLING MATERIAL
ENDOSEAL[®]
Bioceramic Sealer




Proper Setting time

2 - 4 hours



Biocompatibility

Promotes tissue healing without irritating periodontal tissue, thanks to 3 bioactive core ingredients.



Dimensionally Stable

Achieves stability with the lowest solubility compared to market competitors.



Low Film Thickness

Ensures a gap-free, consistent, and effective seal within the root canal.



Superior Radiopacity

7.5mmAl




Aesthetics

Against discoloration



Promotes Bone Generation and Rapid Healing


By Maruchi's patented 3 core formulation.



Initial



1year F/U



Initial



5 mon F/U



Retreatability

Optimized compressive strength allows easy removal using standard hand or NiTi files.

Maruchi's Endo Treatment Unique Solution

Bioceramic Seal

With Optimized Root Repair and Filling



Disinfection

With Effective Medicament and Ultrasonic Irrigation Device



NiTi Rotary Files

Featuring Next-Generation Heat Treatment Technology



Enhance Your Workflow
With Maruchi's Advanced Endo Solutions

We just mimic nature



MARUCHI
www.endoland.com

Tel +82-33-734-0330
Fax +82-33-746-2804
info@endocem.com

101, Donghwagongdan-ro, Munmak-eup, Wonju-si,
Gangwon-do, 26365, Republic of Korea
강원도 원주시 문막읍 동화공단로 101